

5

Date 14/8/19

Expt. No. 15.

Page No. 33




## Volumetric analysis

Aim : (i) To prepare 250 ml of  $M_{10}$  Mohr's salt solution (ferrous ammon. sulphate)  
Using this sol<sup>n</sup> find out the molarity and strength of given  $KMnO_4$  solution.




Apparatus Required : Watch glass, 250 ml conical flask, pipette funnel, weighting machine, titration flask, burette stand.

Chemicals Required : Mohr salt, water, conc  $H_2SO_4$ , dil.  $H_2SO_4$ ,  $KMnO_4$ .

Theory : Redox titration  
Molecular equations



Ionic equations



Teacher's Signature \_\_\_\_\_

Burette : 50 ml =  $\text{KMnO}_4$  = oxidising agent

Pipette : 10 ml = Mohr salt = reducing agent

End point : From colorless to very light pink permanent color

Indicator :  $\text{KMnO}_4$  = self-indicator

### Procedure

Preparation of 250 ml of M/10 Mohr salt solution

Molarity =  $\frac{\text{no. of moles of Mohr salt}}{\text{Vol. of solution in litres}}$

$$M = \frac{\text{mass of Mohr salt}}{\text{molar mass of Mohr salt}} \times \frac{1000}{V_{\text{soln}} (\text{ml})}$$

$$\frac{1}{10} = \frac{w}{392} \times \frac{1000}{250}$$

$$w = 9.8 \text{ g}$$

Weight of empty watch glass = 22 g

Weight of watch glass + Mohr salt = 31.8 g

Weight of Mohr salt = 9.8 g

- 1) Transfer the weighed Mohr salt from watch glass into a clean 250 ml measuring flask using a funnel.
- 2) Wash the watch glass and funnel with distilled water to transfer the particles sticking to it into the measuring flask and add 2-3 ml of conc  $\text{H}_2\text{SO}_4$  to it.

Teacher's Signature \_\_\_\_\_

- 3) Swirl the measuring flask till solid Mohr salt dissolves.
- 4) Add enough water just upto 250 ml.
- 5) Clamp a rinsed burette on burette stand.
- 6) Pipette 10 ml of Mohr salt sol<sup>n</sup> in a titration flask. Add add 1 full test tube (10 ml) of dil.  $H_2SO_4$  to it.
- 7) Place the titration flask under the nozzle of burette.
- 8) Note the initial reading of burette and run out the solution.
- 9) Stop addition of the sol<sup>n</sup> when the end point is reached. Note the reading.
- 10) Repeat the steps till three concordant readings.

### Observation Table

S.No	Initial reading of burette	Final reading of burette	Volume of $KMnO_4$ used.
1.	0	10	10 ml
2.	10	20.1	10.1 ml
3.	20.1	30.1	10 ml
4.	30.1	40.1	10 ml

CONCORDANT VOLUME = 10 ml

### Calculation

formula used ;  $a_1 M_1 V_1 = a_2 M_2 V_2$

(Mohr salt)                      ( $KMnO_4$ )

$$1 \times \frac{1}{10} \times 10 = 5 \times M_2 \times 10$$

$$M_2 = 0.02 \text{ mol/L}$$

Teacher's Signature \_\_\_\_\_



$$\begin{aligned}\text{moles of } \text{KMnO}_4 \text{ sol}^n &= 0.02 \text{ mol/L} \\ \text{strength of } \text{KMnO}_4 \text{ sol}^n &= M \times \text{molar mass of } \text{KMnO}_4 \\ &= 0.02 \times 158 \\ &= \underline{3.16 \text{ g/L}}\end{aligned}$$

### Result

molarity of given  $\text{KMnO}_4$  is 0.02 mol/L  
strength of given  $\text{KMnO}_4$  is 3.16 g/L

### Precautions

- 1) Wash the apparatus with water, and rinse it with the given solution.
- 2) Burette scale reading and eye should be at same level.
- 3) Burette should be clamped vertically in the burette stand.

*Pre L.*  
*26/8*

Next week  
keep it up.

Teacher's Signature \_\_\_\_\_